TTBIO *Bettering your life!*

Automatic Handpiece Maintenance Unit

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TTcare 100+ Instruction Manual

Thank you for purchasing the TTcae 100+, a dental handpiece maintenance unit.
The TTcae 100+ lubricates and internally cleans dental handpieces easily and effectively.
For optimum safety and performance, read this manual thoroughly before using the unit and pay close attention to the warnings and notes.
Keep this manual in a handy place for future reference.

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4. Setup 08

5. Operation 10

6. Maintenance and Replacement Parts 13

7. Maintenance and Inspection 16

8. Trouble shooting and Error Code 17

9. Tools and Consumables 19

10. Service Contacts 20

11. Information on Electromagnetic Compatibility 21

PRECAUTIONS

Most operation and maintenance problems result from insufficient attention being paid to basic safety precautions and not being able to foresee the possibilities of accidents. Problems and accidents are best avoided by foreseeing the possibility of danger and operating the unit in accordance with the manufacturer's recommendations. First thoroughly read all precautions and instructions pertaining to safety and accident prevention; then, operate the equipment with the utmost caution to prevent either damaging the equipment itself or causing injury.

Note the meaning of the following symbols and expressions:

	WARNINGS	This warns the user of danger of death, serious injury or total equipment damage and failure or fire.
	CAUTIONS	This identifies methods not to be used or purposes which the instrument is not suited for.
	NOTES	This alerts the user to the risk of light to medium injury or equipment damage.

■ **Note**
 This alerts the user of important points concerning operation.
 The user (i.e. clinic, hospital etc.) is responsible for the management, maintenance and use of medical devices. Also this equipment must not be used by anyone except legally qualified dentist or doctors.

Do not use this equipment for anything other than its specified purpose.
--

■ **Caution:**
 Federal law restricts this device to sale by or on the order of a dentist (for U.S.A.).

ATTENTION

TTBIO CORP. will not be responsible for accidents, instrument damage, or

- (1) repairs made by personnel not authorized by TTBIO CORP.
- (2) any changes, modifications, or alterations of its products.
- (3) maintenance or repairs using parts or components other than those specified by TTBIO CORP. and other than in their original condition
- (4) operating the instrument in ways other than the operating procedures described in this manual or resulting from the safety precautions and warnings in this manual not being observed
- (5) workplace conditions and environment or installation conditions which do not conform to those stated in this manual such as improper electrical power supply
- (6) fires, earthquakes, floods, lightning, natural disasters, or acts of God.

1. Warnings, Cautions and Notes

WARNINGS

- Explosion Hazard. Do not use near open flames or other ignition sources.
- Health Hazard. The area must be well ventilated. Breathing fumes could damage your health.
- The equipment must be properly grounded.
- Make sure the ground is properly connected.
- To prevent the spread of grave, life-threatening diseases like HIV and hepatitis B, autoclave handpieces after performing regular maintenance using the TTcare 100+.
- Avoid the risks of electrical shock, equipment damage and fire during an electrical storm: Turn the TTcare 100+ off and do not touch it or its cord.
- Wear surgical gloves, mask and protective eyewear to operate and clean the TTcare 100+.

CAUTIONS

Electromagnetic wave interference caused by cellular phones, transceivers, remote controls and similar transmission devices could cause the equipment to operate randomly. All devices which transmit electromagnetic waves located near the work area should be turned off.

NOTES

- Set the TTcare 100+ on a level and stable surface.
- Depending on the type of handpiece and coupling, there could be a release of oil mist. In this case, lower the pressure to 4 bar.
- Oil will leak or spray out if the cans are not properly installed.
- Make sure the head is flat against the nozzle guide. If it is tilted, the nozzle might be bent or oil might spray out.
- Oil will spray out of the chuck nozzle if the handpiece is removed too soon (before the beep).

- For some handpieces, an oil mist may be released if a burr is not inserted.
- Do not damage the o-ring on the coupling when attaching a handpiece. This could result in oil spraying out or poor performance.
- Do not select a coupling that does not have a handpiece attached to it. This would result in oil spraying out when the Start Switch is pressed.
- If the unit will not be used for a while, unplug it and close the main air valve.
- Please insure you use the correct adaptor for the instrument.
- Do not use couplings made for maintenance equipment other than the TTcare 100+. This could result in release of oil mist or poor handpiece maintenance.
- Dispose of used oil pads according to regulations for medical waste.
- Do not fail to push the tray all the way in. Otherwise, oil could leak out or oil mist might be released.

2. Specifications

Name	TTcare 100+
Rating	AC 100-240 V 50/60 Hz
Power Consumption	60VA max.
Fuse	250V 2A Slow Blow Type
Air Pressure	4~6 bar
Recommended Air Pressure	4 bar
Air Flow Rate	40 NL/ min
Weight	Approx. 9 kg
Size (including regulator)	Width 295 × Height 385 × Depth 295 mm
Protection category (Dust-proof; water-proof)	IP20

Operating environment

Ambient temperature	Permitted in interior rooms
Temperature	15 to 40°C (59 to 104°F)
Humidity	25 to 90 %
Max. altitude	2000 m (6560 feet)

Transportation and storage conditions

Temperature	-20 to 70°C (-4 to 158°F)
Humidity	5 to 95 %
Atmospheric pressure	700 to 1060 hPa

Product Description

Used to maintain optimum performance and prolong working life of dental handpieces.
 Delivers oil and air automatically to handpiece.
 Used after dental treatment and before autoclaving.

Disposal

The wastes incurred are to be recycled or disposed of in a way that is harmless for human beings and the environment; in doing so, the national valid regulations are to be observed.
 Disposal of equipment and accessories at the end of their service lives:
 On the basis of EC Directive 2002/96/EC on Waste Electrical and Electronic Equipment, we would like to point out that this product is currently in compliance with the labelling requirements but is not yet subject to the disposal requirements of this Directive. However, the unit may be disposed of in Europe in special waste management centers. Additional information can be obtained from the manufacturer or your dental supplier.

Meaning of the Symbols

• Device Label

TTcare 100+
Automatic Handpiece Maintenance Unit

TTBIO CORP.
2F, No.7, 6th Road Industry Park Taichung, Taiwan 40755

Product code : 08090012XXX
SN XXXXXXXXXXXX

INPUT :AC 100~240V 50/60 Hz
60VA max.

AIR PRESSURE: 4~6 bar
FUSE: 250V 2A Slow Blow Type

WARNING

- The blister can be cleaned with 30°C±5°C tap water or 60% to 70% alcohol, non-high-temperature cleaning.
- The blister can be used repeatedly, but should be changed once broken.

FLAMMABLE

- Spray cans will explode if they get too hot.
 - Do not expose the spray cans to heat.
 - Only use approved spray can.
- * PREVENT EXPLOSION AND IGNITION OF SPRAY CANS!**

• Appearing on the device

	Observe operator's manual	SN	Serial Number
	Date of manufacture		Warning
	Disposal information, see important information " Disposal"		Flammable

IP20	Protection category Dust-proof; water-proof		Eye protection must be worn
-------------	---	--	------------------------------------

• Carton Label

TTcare 100+
Automatic Handpiece Maintenance Unit

TTBIO CORP.
2F, No.7, 6th Road Industry Park
Taichung, Taiwan 40755

Product Code:08090012 SN XXXXXXXXXXXX

Storage and transport requirements

Temperature range

Air Pressure

Humidity

Transport upright

Fragile

Stacking restriction

Keep dry

Waste Disposal

• Appearing on the packaging

	Fragile		Air pressure
	Keep dry		Humidity
	Transport upright		Quantity
	Stacking restrictions		Manufacturer
	Temperature range		



Main Switch

	OFF		ON
--	------------	--	-----------


05

3. Parts Identification

■ Introduction

■ Components

<p>Main Unit</p> 	<p>Blister (CODE No. 0PL0350)</p> 	<p>Absorption pads (CODE No. 0PL0351)</p> 
<p>Oil Tray (CODE No. 0PL0359)</p> 	<p>Air Tube (CODE No. 0PL0360)</p> 	

* TTcare 100+ contains above standard components.
* Optional components see chapter 9.

4. Setup

■ Note


Leave at least 5 cm (2 in) of free and open space around the TTcare 100+.

(1) Air Connection

1. Connect the air tube to the air connection and the main air supply.
Turn on the valve for the main air supply.

■ Note


The main air supply must be clean and oil-free.
Make sure the tube is securely connected.

2. Set the air regulator dial between 4 and 6 bar.
Recommended Air Pressure is 4 bar.


■ Note

Less than 4 bar or higher than 6 bar will result in an error code.

(2) Power Cord

- Connect one end of the power cord to the main unit and plug the other end into the electrical supply.


(3) Adapter

- Fit the adapter into the nozzle and fasten the nut.

■ Note

Make sure the nut is securely tightened. Otherwise oil and air will escape and cause incorrect maintenance.


(4) Tray

- Place two Absorption pads in the tray and slide it into place.

■ Note

Slide the tray all the way in; otherwise, it might leak.

(5) Spray Cans


1. Open the side door.
Press the PUSH sign indicated in the photo to open the door.

■ Note

Do not use Acidity Cleaner which is used for air bearing handpieces only. It could damage both the TTcare 100+ and the handpiece.


2. Ensure the lock lever is on the left, place the spray can on the stand shown in the photo.

3. Shift the lock lever to the right to secure the spray can.

■ Note



Ensure the nozzle is properly connected.
Make sure the lock lever is all the way to the right, otherwise the oil will not come out.
Gently pull the spray can to ensure it is secured.

! NOTE

Oil will leak or disperse if the cans are not properly installed.

4. Gently lift the spray can and ensure the nozzle is secured.

5. Use a hex wrench to adjust the screw on the spray can stand, if necessary.

■ **Note**


This can make up for a difference in height of 5mm.

NOTE

- It is recommended to use TT BIO or TTlube spray can for the TTcare 100+.
- The spray can nozzle and spray can stand are designed to offer compatibility with other spray can brands but do not mean all applicable.

6. Close the side door.

5. Operation


■ **Panel display:**

■ **Operating instruction:**

Icons	Functions	Illustration	
	Power On/Off		OFF
			ON
	Selection of handpiece 1	Handpiece 1 is selected when the lamp lights up.	
	Selection of handpiece 2	Handpiece 2 is selected when the lamp lights up.	
	Selection of handpiece 3	Handpiece 3 is selected when the lamp lights up.	
	Start / Set / Chuck maintenance switch	1. Handpiece maintenance when door is closed 2. Chuck maintenance when door is open 3. Function setting by long press and wait for the beep	
	Spray can change	Replace spray can when lamp lights up	

■ **Operating instruction:**

■ Operation maintenance:

■ Note

Don't insert burr for processing maintenance!

(1) Chuck maintenance:

Turn on the power switch. Open the front door.

will show on LCD Display during process.

Place the chuck side of head to the chuck nozzle.

Press to begin cleaning and lubricating the chuck.

A beep sound meaning the procedure is completed.

(2) Handpiece maintenance:

Turn on the power switch.

Open the front door.

Connect the handpiece to its coupling with the chuck facing back.

Make sure the handpieces are connected to the adapter until it clicks into place.

Close the front door.

will show on LCD Display during process.

Press / / to select handpieces which are connected.

Press to lubricate the handpieces in sequence order.

A beep sound meaning the procedure is completed.

■ Note

Pressing any button or opening the front door will immediately stop the maintenance procedure.

■ Functional setting:

Select functional setting mode by pressing for approx 2 sec.

- (1) will show on LCD Display meaning compressed air is on.

Turn compressed air On/Off by pressing .

- (2) will show on LCD Display by pressing shortly.

Maintenance time setting:

Select 15 sec by pressing / Select 20 sec by pressing / Select 30 sec by pressing .

- (3) will show on LCD Display meaning buzzer sound is on by pressing shortly.

Select buzzer sound On/Off by pressing .

- (4) will show on LCD Display meaning auxiliary motor is on by pressing

shortly. Select Auxiliary motor On/Off by pressing .

- (5) will show on LCD Display by pressing shortly. Press / / to select handpieces which are connected to run optional auxiliary rotation.
- (6) will show on LCD Display by pressing shortly.

■ Compressed air setting:

Adjust air pressure by turning between 4.0-6.0 bars. Recommend to set air pressure at 4.0bars for TTcare 100+. Save and leave functional setting mode by pressing shortly.

■ Oil amount adjustment :


- (1) Indicator lights and are on when the Lubrication Unit is turned on. Press for three seconds until it enters Oil amount adjustment mode, LCD will display:

- (2) There are 11 adjustable oil amount from PU-00 to PU-10. Press for increasing or for decreasing. PU-00 is the default setting, PU-10 has the largest oil amount.

Choose the PU value required, press button to complete setup. Lubrication unit will automatically reboot after setting adjustments.

6. Maintenance and Replacement Parts

■ Replace Fuse

- (1) Press fuse holder by slotted wrench to loose the holder.
- (2) Remove fuse holder.
- (3) Replace new fuse (250V2A Slow-Blow type) into fuse holder
- (4) Insert fuse holder.

■ Remove and Clean Blister


⚠ WARNING


- The blister can be cleaned with 30°C±5°C tap water or 60% to 70% alcohol, non-high-temperature cleaning.
- The blister can be used repeatedly, but should be changed once broken.

■ Remove and Clean Front Door:

- (1) Push left pivot upward and inward.
- (2) Push right pivot upward and inward.

- (3) Pull front pin in the direction of arrow shown in the photo.
- (4) Remove front door.

The front door can be cleaned with 30 ° ± 5 ° C tap water or 60% to 70% alcohol, non-high-temperature cleaning.

■ Cleaning the Tray:

- (1) Take the tray out in the direction of arrow shown in the photo.
- (2) Clean and replace new absorption pad.

⚠ CAUTION

- The tray cannot be autoclaved.
- Autoclaving will cause it to rust.
- Do not use disinfectants that contain chlorine-based cleaners or benzalkonium chloride or chlorinated aromatics.
- If this tray is not cleaned, oil will overflow from the tray.

! NOTE

- Dispose of used absorption pad according to regulations for medical waste.
- Push the tray all the way in, otherwise, oil may leak and oil mist might be dispersed.

■ Cleaning Couplings :

Loosen the nut and take the coupling off.

■ Note


- (1) Remove it carefully without damaging the o-ring.
- (2) Couplings cannot be autoclaved.

Clean it with 60% to 70% alcohol.

An o-ring could slip out of its groove. Make sure all o-rings are properly seated in their grooves before attaching a handpiece.

■ Note

Do not connect handpiece to the coupling when o-ring is not properly fitted in its groove. Possible damages to the o-ring and stuck in the coupling may make it difficult for disassembling handpiece.

7. Maintenance and Inspection

* Send the product every 2 years for a service check.

The safety checks in different countries can vary in compliance with country specific regulations and requirements for medical devices. The national valid regulations are to be observed.

* Repair and maintenance work- part from the activities described in these operating instructions-may be performed only by qualified technical staff.

* Expected service life 5 years.

• In the event of modification by third parties, existing medical device licences become null and void.

• Only use original parts and spare parts.

• The user(i.e. clinic, hospital etc.) is responsible for the management, maintenance and use of medical devices.

■ Regular Inspection

* Maintenance and inspection are generally considered to be the duty and obligation of the user, but if, for some reason, the user is unable to carry out these duties, they may rely on a qualified medical device serviceman. Contact your local dealer or the TTIBO CORP. for details.

(1) Power Supply Cord

Inspect visually for wear and broken wires.

(2) Main and Operation Switches

Turn on main switch and ensure main lamp lights up.

Check that the unit operates correctly.

(3) Absorption pads

Replace absorption pads(when the tray is filled by oil)..

• Regularly Required Replacement Parts

Fuse

• Other Replacement Parts

Absorption pads, Blister

8. Trouble Shooting and Error Code

■ Trouble Shooting

Breakdown	Reason	Solution
Main power lamp does not light up	Power supply cord not properly plugged into supply socket	Plug in properly
	Power supply cord not properly plugged into TTcare100+	
Does not start when Start Switch is pressed. (1) Double beep sound (2) Spray can Empty Lamp is blinking	Front door may be open	Close front door
	Empty can	Replace can
	Can is not properly installed	Properly install the can
	Compressor air is not connected	Connect air
Unusual odor. Surfaces are unusually dirty.	Air regulator is not set properly	Set regulator properly
	Use without connecting handpiece	Select only couplings that have a handpiece connected
	Absorption pads is saturated	Clean TTcare100+ and replace filter
	Can is not properly installed	Properly install the can
Cannot install spray can. Oil leaks from spray can.	Air regulator is not set properly	Set regulator properly
	Bent spray can nozzle pin	Replace can
Oil mist is released. Equipment surfaces get oily	Wrong stand for spray can	Use proper stand for can
	Type of handpiece or coupling	Depending on the type of handpiece and coupling, there could be a release of oil mist. In this case, lower the pressure to 4 bar
	Air pressure too high	Reduce air pressure, but not less than 4 bar
	No burr inserted in straight attachment.	Refer to the user instructions for the handpiece. For some handpieces, an oil mist may be released if a burr is not inserted

■ Error Code

Error Code	Description	Reason	Exclusion
Error 01	Keyboard Error	Function keypad Error	Contact your dealer
Error 02	Oil supply motor Error	Faulty Motor or internal circuit	Contact your dealer
Error 03	Auxiliary rotation motor Error	Faulty Motor or internal circuit	Contact your dealer
Error 04	Pressure input too low	Air pressure lower than preset	Increase air pressure renew air setting(refer to chapter 4)
Error 05	Pressure input too high	Air pressure higher than preset	Decrease air pressure renew air setting(refer to chapter 4)
Error 06	Function Memory Access Error	Faulty Main IC	Contact your dealer
Error 07	Oil empty	Oil empty	Replace a new can of oil or Contact your dealer
Error 08	Electromagnetic valve Error	Faulty Electromagnetic valve	Contact your dealer
Error 09	Circuit is over temperature	Faulty Motor or internal circuit	Contact your dealer
Error 10	Circuit detect Error	Faulty Main IC	Contact your dealer

9. Tools and consumables

Presentation	Material summary	Part No.
	KaVo service adapter	0PL0345
	W&H service adapter	0PL0346
	NSK service adapter	0PL0344
	Bien-Air service adapter	0PL0348
	SIRONA service adapter	0PL0347
	E-Type adapter	0PL0349
	Spray can	08090011
	Blister	0PL0350
	Absorption pads	0PL0351

10. Service Contacts

■ For repair or other types of service contact your local dealer or TT BIO CORP.

■ Disposal of Medical Devices

Any medical devices which could possibly be contaminated must be first decontaminated by the responsible doctor or medical institution and then be disposed by an agent licensed and qualified to handle medical and industrial waste.

11. Information on Electromagnetic Compatibility

Guidance and manufacturer's declaration-electromagnetic emissions

The TTcare 100+ is intended for use in the electromagnetic environment specified below. The customer or the user of the TTcare 100+ should assure that it is used in such an environment.

Emission Test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The TTcare 100+ uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Group B	The TTcare 100+ is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Group A	
Voltage fluctuations /flicker emissions IEC 61000-3-3	Compliance	

Guidance and manufacturer's declaration-electromagnetic immunity			
The TTcare 100+ is intended for use in the electromagnetic environment specified below. The customer or the user of the TTcare 100+ should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge(ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or eramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC 61000-4-4	± 2kV for power supply lines ± 1kV for input/output lines	± 2kV for power supply lines Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1kV line(s) to line(s) ± 2kV line(s) to earth	± 1kV differential mode Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT(>95% dip in UT) for 0,5 cycle 40% UT(60% dip in UT) for 5 cycles 70% UT(30% dip in UT) for 25 cycles <5% UT(>95% dip in UT) for 5 s	<5% UT(>95% dip in UT) for 0,5 cycle 40% UT(60% dip in UT) for 5 cycles 70% UT(30% dip in UT) for 25 cycles <5% UT(>95% dip in UT) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the TTcare 100+ requires continued operation during power mains interruptions, it is recommended that the TTcare 100+ be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	The TTcare 100+ power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
<p>■ Note: UT is the a.c. mains voltage prior to application of the test level.</p>			

Guidance and manufacturer's declaration-electromagnetic immunity			
The TTcare 100+ is intended for use in the electromagnetic environment specified below. The customer or the user of the TTcare 100+ should assure that is used in such and environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 KHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the TTcare 100+ including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance: $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P}$ 80MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800MHz to 2,5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3 V/m 80MHz to 2,5 GHz	3 V/m	
<p>■ Note1: At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>■ Note2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>■ a: Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the TTcare 100+ is used exceeds the applicable RF compliance level above, the TTcare 100+ should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the TTcare 100+.</p> <p>■ b: Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Recommended separation distance between portable and mobile RF communications equipment and the TTcare 100+			
The TTcare 100+ is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the TTcare 100+ can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the TTcare 100+ as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter M		
	150 kHz to 80 MHz d = 1,2√P	80 MHz to 800 MHz d = 1,2√P	800 MHz to 2,5 GHz d = 2,3√P
0.01	0.02	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
<p>■ Note1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.</p> <p>■ Note2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

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Recommended separation distance between portable and mobile RF communications equipment and the TTcare 100+			
The TTcare 100+ is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the TTcare 100+ can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the TTcare 100+ as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter M		
	150 kHz to 80 MHz d = 1,2√P	80 MHz to 800 MHz d = 1,2√P	800 MHz to 2,5 GHz d = 2,3√P
0.01	0.02	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

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■ **Note1:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

■ **Note2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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The TTcare 100+ is intended for use in the electromagnetic environment specified below.
The customer or the user of the TTcare 100+ should assure that it is used in such an environment.

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Guidance and manufacturer's declaration-electromagnetic immunity			
The TTcare 100+ is intended for use in the electromagnetic environment specified below. The customer or the user of the TTcare 100+ should assure that is used in such and environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 KHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the TTcare 100+ including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance: $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P}$ 80MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800MHz to 2,5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3 V/m 80MHz to 2,5 GHz	3 V/m	
<p>■ Note1: At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>■ Note2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> <p>■ a: Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the TTcare 100+ is used exceeds the applicable RF compliance level above, the TTcare 100+ should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the TTcare 100+.</p> <p>■ b: Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

**Recommended separation distance between
portable and mobile RF communications equipment and the TTcare 100+**

The TTcare 100+ is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the TTcare 100+ can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the TTcare 100+ as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter M		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3\sqrt{P}$
0.01	0.02	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

- **Note1:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
- **Note2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.