www.ttbio.com

## 1. Product Description

### Tiger 101 Design

- Connection (Borden 2-hole, Midwest 4-hole)
- Head size (T: torque, M: mini)
- Spray hole (Single, Triple)
- Press chucking

## 2. Specifications

	Fixed Connection	
	Tiger 101-T2 / M2	Tiger 101-T4 / M4 / 3T4 / 3M4
Connection	Borden 2-hole	Midwest 4-hole
Spray Hole	Single	Single / Triple
Rotation speed (rpm)	T≥300,000	M≥350,000
Bur size-diameter/length (mm)	Ø1.59~Ø1.60 / T(19-25) ; M(16-21)	
Drive air pressure (bar)	2.2~2.8	
Suggested drive air pressure (bar)	2.5	
Water flow rate (ml/min)	>50	
Water pressure (bar)	0.8~2.0	
Air pressure (bar)	1.0~3.0	

- ※ The ISO 2/4 holes is defined in accordance with ISO 9168:2009
- ※ Air pressure must be higher than water pressure. ※ T: torque head M: mini head
- ※ Pressure values are measured by multi gauge.

Air compressor and other devices used by the handpiece must be in accordance with ISO requirement. (MUST set air pressure higher than water pressure before operation.)

- Using contaminated water and compressed air containing moisture, impurity or oil to drive air could shorten the lifetime of the cartridge/rotor, and block the spray nozzle.
- If NOT follow the required configuration, the risks of damaging the devices and injury to people could occur.
- Handpiece is reusable. After unpacking from package or before operating, please proceed the sterilization from step 7 to 9 by following the procedures: cleaning>lubrication>sterilization.

## 3. Safety Notes

- Once handpiece is damage (ie deformed caused by hit) or improperly use (ie sterilized above suggested temperature), be sure not to continue using it.
- Before operation, always make sure the product and its accessories meet the specified requirements (ie drive air pressure, lubrication...etc).
- Please contact our dealer or us with proof of purchase, whenever you have problems such as abnormal vibration, exceeding temperature or unusual noise...etc.
- Avoid contact between the instrument head and soft tissue (risk of burning due to the press button heating up).

## 4. Installation


Fig.1

Fig.2



Fig.3

To install handpiece into ISO 2/4-hole hose(\*), please follow the steps illustrated from Fig.1 to Fig.3. Make sure the handpiece is connected to the hose tightly. If there is any leakage of water or air, please repeat this installation procedure.

\* The holes hose is defined in accordance with ISO 9168 : 2009

## 5. Changing Bur

- Follow arrow direction shown in Fig. to take off or install the bur.
- Put on protective gloves to reduce the risk of interactive infection, especially not to touch the used burs without proper protection.
- Make sure bur is tightly secured before operation, otherwise risks of injury to doctors, patients, or the third party could occur.
- Bur must be complied with ISO requirement as described in the Spec. Burs not in good conditions such as being bent and damaged, or are deviated from ISO standards and in improper length may decrease the lifetime of the handpiece or even cause severe consequence such as injury to people.
- Do not press the button while bur is rotating, otherwise it could harm internal parts.


## 6. Replacement of Rotor

- Please contact local dealers and let the technicians replace the rotor.
- Do not use unsuitable rotors or spare parts, otherwise could damage the handpiece and result in severe consequence during treatment.

## 7. Cleaning

- Put on protective gloves to detach the handpiece from the hose or coupling.
- Use clean cloth with mild cleaning solution to clean the handpiece surface carefully, especially handpiece groove, and remove the residues.
- Do not immerse handpiece or cartridges in disinfectant solution or water; Not to be cleaned by ultrasonic process.
- Cleaning solution must be non-erosive and harmless to human. Solution with 80% ethyl alcohol is recommended. Please also follow recommendation of the cleaning solution manufacturer for suggested usage.
- Clean the spray nozzle with the nozzle cleaner needle when it is blocked to prevent burn caused by overheat (see above figure).

## 8. Lubrication and Maintenance

- Use dental handpiece lubrication oil for maintenance, any unspecified oil not intended for lubricating handpieces could shorten the lifetime and harm human body.
- To lubricate 101 series handpiece, put Spindle Nozzle into nozzle of the lubrication oil can and inject the oil into drive air tube for 2~3 seconds (Fig.4).
- Be sure to lubricate the handpieces prior to sterilization and run the handpieces for 20~30 seconds to emit residual oil and dirt.
- To lubricate handpiece chucks, put Spindle Nozzle into nozzle of the lubrication oil can and inject lubrication oil into the chuck (Fig.5) for 1 second. It's recommended to lubricate chuck at least once a week and then proceeds step 2.
- After unpacking and prior to first use, be sure to lubricate with oil and run the handpieces in low rotation speed for 30 seconds and in full speed for 20~30 seconds to emit residual oil.


Fig.4

Fig.5

## 9. Sterilization

TTBIO suggested sterilization in accordance with EN13060

- Place handpiece in single-use sterile bag and proceed steam prevacuum sterilization. Do not open the sterile bag before use in order to avoid contact infection.
- Use qualified sterilization equipment specified by the authority/government for sterilization procedures.
- Do not use chemical steam nor immerse handpiece or cartridge/rotor in disinfectant solution, otherwise could damage the handpiece and cause severe consequences.
- Sterilization holding time a minimum of 4 minutes at 134°C. (drying time is 15 minutes)

## 10. Disclaimer

- Information in this operation manual were carefully reviewed. We cannot assume any liability for any damage arising from improper use of the handpiece or from use of this manual.
- Products are subject to change without notice.
- This product is only for dental handpiece treatment and use for other purpose may damage the handpiece and cause severe consequence.
- Please refer to the bar code on the packing box for manufacturing and expiration date.

CAUTION: U.S. Federal law restricts the use of this device to licensed professionals. (Rx Only)

## 11. Quality and Service

All of our products have passed strict QC procedures before delivery. Please note that improper use or not using the recommended parts and accessories or modifying our products without our approval may shorten the lifetime of the product or cause damage. Please contact us or dealers with proof of purchase for consultation whenever you have any questions.

### Dealer

Manufactured by  
TTBIO CORP.  
2F.,NO.7, 6TH ROAD INDUSTRY PARK  
TAICHUNG, TAIWAN 40755  
TEL:886-4-23595958  
FAX:886-4-23596893  
E-mail:medical@ttbio.com  
Website:www.ttbio.com

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