CE
0086

www.ttbio.com

1. Product Description

Tiger 500⁺ Design

- Omni-directional lighting system
- Titanium coating
- Quadruple spray hole
- Push button
- Grip profile

2. Specifications

Model	Tiger 500	Tiger 500K	Tiger 500N
Connection	TTBIO® EZclick	TTBIO® EZclick-KL/ KaVo® Multiflex®	TTBIO® EZclick-NL/ NSK® PLT-CL®
Head size & light	M ⁺ / ML ⁺ / T ⁺ / TL ⁺ (*)		
Rotation speed (rpm)	T≥300,000 M≥350,000		
Bur size-diameter/length (mm)	Ø1.59~Ø1.60 / T(19~25) ; M(16~21)		
Drive air pressure (Bar)	2.6~3.0		
Recommended Drive air pressure (Bar)	2.8 (***)		
Spray Volume (ml/min)	> 50		
Water pressure (Bar)	0.8~2.0(**/***)		
Chip air pressure (Bar)	1.0~3.0		

* M= Mini Head; T= Torque Head; L= Light Support

** Chip Air Pressure must be greater than Water Pressure and both values need to be set at the same time

*** Measured by multi-gauge

3. Safety Notes

1. △ Once handpiece is damaged (ie deformed caused by hit) or improperly used (ie sterilized above suggested temperature), be sure not to continue using it.
2. Before operation, always make sure the product and its accessories meet the specified requirements (ie drive air pressure, lubrication...etc).
3. Please contact our dealer or us with proof of purchase, whenever you have problems such as abnormal vibration, exceeding temperature or unusual noise...etc.
4. Avoid contact between the instrument head and soft tissue (risk of burning due to the press button heating up).

4. Installation

Fig. 1 Install EZclick connector

Fig. 2 Rotate the connector to install


Fig. 3 Strip EZclick connector

To install or strip handpiece, please follow the steps illustrated from Fig. 1 to Fig. 3 (Please refer to the coupling manufacturers' instructions.)

5. Changing Bur

1. Follow the arrowhead direction shown in Fig. to take off or install the bur.
2. Put on protective gloves to reduce the risk of interactive infection, especially not to touch the used burs without proper protection.
3. Make sure bur is tightly secured before operation, otherwise risks of injury to doctors, patients or the third party could occur.
4. Bur must be complied with ISO requirement. Burs not in good conditions such as being bent and damaged, or are deviated from ISO standards and in improper length may decrease the lifetime of the handpiece or even cause severe consequences such as injury to people.
5. Do not press the button while bur is rotating, otherwise it could harm internal parts.

6. Replacement of Rotor

1. Please contact local dealers and let the technicians replace the rotor.
2. Do not use unsuitable rotors or spare parts, otherwise could damage the handpiece and result in severe consequences during treatment.

7. Cleaning

1. Put on protective gloves to detach the handpiece from the coupling.
2. Use clean cloth with mild cleaning solution to clean the handpiece surface carefully, especially handpiece groove, and remove the residues.
3. △ Do not immerse handpiece or cartridges/rotor in disinfectant solution or water, not to be cleaned by ultrasonic process.
4. △ Cleaning solution must be non-erosive and harmless to human. Solution with 80% ethyl alcohol is recommended. Please also follow recommendation of the cleaning solution manufacturer for suggested usage.
5. Clean the spray nozzle with the nozzle cleaner needle when it is blocked (see above figure)

8. Lubrication and Maintenance

1. △ Use dental handpiece lubrication oil for maintenance, any unspecified oil not intended for lubricating handpieces could shorten the lifetime and harm human body.
2. To lubricate 500⁺ series handpiece, put Spray Nozzle into nozzle of the lubrication oil can (Fig. 4). Connect handpiece with the nozzle (Fig. 5) and inject lubrication oil for 2~3 seconds.
3. Be sure to lubricate the handpieces prior to sterilization and run the handpieces for 20~30 seconds to emit residual oil and dirt.
4. To lubricate handpiece chucks, put Shaft Spray Nozzle into nozzle of the lubrication oil can and inject lubrication oil into the chuck (Fig. 6) for 1 second. It's recommended to lubricate chuck at least once a week and then proceed step 2.
5. After unpacking and prior to first use, be sure to lubricate with oil and run the handpieces in low rotation speed for 30 seconds and in full speed for 20~30 seconds to emit residual oil.


Fig. 4


Fig. 5


Fig. 6

9. Sterilization

TTBIO suggested sterilization in accordance with EN13060

1. △ Place handpiece in single-use sterile bag and proceed steam prevacuum sterilization. Do not open the sterile bag before use in order to avoid contact infection.
2. △ Use qualified sterilization equipment specified by the authority/government for sterilization procedures.
3. △ Do not use chemical steam nor immerse handpiece or cartridge/rotor in disinfectant solution, otherwise could damage the handpiece and cause severe consequences.
4. △ Sterilization holding time is a minimum of 4 minutes at 134°C. (drying time is 15 minutes)

10. Disclaimer

1. Information in this operation manual were carefully reviewed. We cannot assume any liability for any damage arising from improper use of the handpiece or from use of this manual.
2. Products are subject to change without notice.
3. This product is only for dental handpiece treatment and use for other purpose may damage the handpiece and cause severe consequences.
4. Please refer to the bar code on the packing box for manufacturing and expiration date.


KaVo® & MULTiflex® are the registered trademarks of Kaltenbach & Voigt GmbH CO., Germany. NSK® & PLT-CL® are the registered trademarks of NAKANISHI INC. Japan.

CAUTION: U.S. Federal law restricts the use of this device to licensed professionals.

11. Quality and Service


All of our products have passed strict QC procedures before delivery. Please note that improper use or not using the recommended parts and accessories or modifying our products without our approval may shorten the lifetime of the product or cause damage. Please contact us or dealers with proof of purchase for consultation whenever you have any questions.

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