
Vigor Air Motor


Operation Manual

1. Product Description

Vigor A65 Air Motor

Vigor A68 Air Motor

2. Specifications

Air motor	Vigor A65	Vigor A68
Hose connection	ISO 9168	ISO 9168
Instrument connection	ISO 3964	ISO 3964
Speed range (rpm)	5,000-20,000 (2.5 Bar) 5,000-22,000 (3.0 Bar)	5,000-20,000 (2.5 Bar) 5,000-22,000 (3.0 Bar)
Spray volume (ml/min)	>50	>50
Water pressure (Bar)	1.5 - 2.5 *	1.5 - 2.5 *
Spray air pressure (Bar)	2.0 - 3.0 *	2.0 - 3.0 *


* Spray Air Pressure must be greater than Water Pressure and both values need to be set at the same time

3. Safety Notes

- (1) The Vigor Air Motor are intended only for use within dentistry as part of dental treatment; any form of misuse or modifications to the product are prohibited and may cause danger.
- (2) Once air motor is damaged (ie deformed cause by hit) or improperly used (ie sterilized above suggested temperature), be sure not to continue using them.
- (3) Before operation, always make sure the product and its accessories meet the specified requirements.
- (4) It is the duty of the user to use only faultless tools, for the correct purpose, and protect patients and third parties from all dangers.
- (5) Higher temperature might cause the sense of oral surface tissue unadaptable. Stop using the handpiece if overheated.


4. Forward / Reverse Rotation of Air Motor

- (1) Forward: Turn the FWD/REV Switch Ring to "F" Mode.
- (2) Reverse: Turn the FWD/REV Switch Ring to "R" Mode.
- (3) When FWD/REV Switch Ring is positioned in the middle between "F" and "R", the air motor stops operating.

5. Water and Chip Air Supply

- (1) Vigor A65 :
 - (A) When water is required at the burr, connect the spray tube attachment to the water supply nozzle of the motor and connect the other end to the water supply nozzle of the Straight or Contra-angle Handpiece.
 - (B) Trim the tube to the required length.
 - (C) The spray tube attachment for Straight handpiece is an option(*0PL0227).

- (2) Vigor A68 :
 - (A) When water is required, connect relevant Straight or Contra-angle Handpiece.
 - (B) When water leak after long use, replace the Oring(*0PL0300).

6. Cleaning

- (1) Use clean cloth with mild cleaning solution to clean the air motor surface carefully, especially groove, and remove the residues.
- (2) Do not immerse air motor in disinfectant solution or water. Do not clean by ultrasonic process.
- (3) Cleaning solution must be non-erosive and harmless to human. Solution with 80% ethyl alcohol is recommended. Please also follow recommendation of the cleaning solution manufacturer for suggested usage.

- (1) Remove the transmission instrument from the air motor and set speed control ring to F or R.
- (2) Remove the air motor from the supply hose.

7. Lubrication

- (1) Put Shaft Spray Nozzle(*0PL0044) into nozzle of the lubrication oil can.
- (2) Inject the oil into the drive air pipe for 2~3 seconds, at the back end of the motor.

- (1) Use lubrication oil to maintain air motor will extend lifetime, any unspecified oil not intended for lubricating air motor could shorten the lifetime and harm human body.
- (2) Be sure to lubricate the air motor prior to sterilization and run the air motor for 20~30 seconds to emit residual oil and dirt.
- (3) After unpacking and prior to first use, be sure to lubricate with oil and run the air motor in low rotation speed for 30 seconds and in full speed for 10~20 seconds to emit residual oil.
- (4) Warnings about lubricant hazards: selected spray can labeling and packaging should comply with other Federal requirements for transport and use.

8. Sterilization

- (1) Place air motor in single-use sterile bag and proceed steam gravity sterilization. Do not open the sterile bag before using in order to avoid contact infection.
- (2) Use qualified sterilization equipment specified by the authority/government for sterilization procedures.
- (3) Do not use chemical steam nor immerse air motor in disinfectant solution, otherwise could damage the handpiece and cause severe consequences.
- (4) Sterilization time is a minimum of 4 minutes at 134°C, and drying time is 15 minutes (according to EN 13060).


9. Disclaimer

- (1) Once this product is damaged (ie deformed caused by hit) or improperly used (ie sterilized above suggested temperature), be sure not to continue using it.
- (2) Before operation, always make sure the air motor and accessories meet operational requirements (ie drive air pressure, water pressure and lubrication, etc)
- (3) Information in this operation manual were carefully reviewed. We cannot assume any liability for any damage arising from improper use of this product or from use of this manual. Please contact our dealer or us with proof of purchase whenever you have problems such as abnormal vibration, exceeding temperature or unusual noise ... etc.
- (4) Products are subject to change without prior notice.
- (5) Products must be operated by qualified practitioners and in dental treatment. Use for other purposes could result in dangers.
- (6) U.S. Federal law restricts the use of this device to licensed professionals. (Rx only)

10. Quality and Service

All of our products have passed strict QC procedures before delivery. Please note that improper use or not using the recommended parts and accessories or modifying our products without our approval may decrease the lifetime of the product or cause damage. Please contact our dealers or us with proof of purchase for consultation whenever you have any questions.

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