
Vigor Series Operation Manual

Vigor S15 Straight Handpiece
Vigor C36 Contra-angle Handpiece
Vigor A65 Air Motor

1. Product

Device description:

The Vigor Series Low speed dental handpieces & Accessories are powered by either low speed air motor or electric micro-motor for teeth cutting, cavity and crown preparation, restorations and polishing teeth.

2. Specifications

Contra-angle	Vigor C36	Air motor	Vigor A65
Transmission ratio	1 : 1	Hose connection	ISO 9168
Motor coupling	ISO 3964	Instrument connection	ISO 3964
Instrument ISO 1797-1 Type 1(ømm)	2.35	Speed range (rpm)	5,000-20,000 (2.5 Bar) 5,000-22,000 (3.0 Bar)
Allowable Max. speed (rpm)	25,000	Spray volume (ml/min)	>50
		Water pressure (Bar)	1.5 - 2.5 *
		Spray air pressure (Bar)	2.0 - 3.0 *

Straight handpiece	Vigor S15
Transmission ratio	1 : 1
Motor coupling	ISO 3964
Instrument ISO 1797-1 Type 2(ømm)	2.35
Allowable Max. speed (rpm)	40,000


* Spray Air Pressure must be greater than Water Pressure and both values need to be set at the same time

3. Safety

- The Vigor series handpieces and air motor are intended only for use within dentistry as part of dental treatment; any form of misuse or modifications to the product are prohibited and may cause danger.
- Once handpiece or air motor is damaged (ie deformed cause by hit) or improperly used (ie sterilized above suggested temperature), be sure not to continue using them.
- Before operation, always make sure the product and its accessories meet the specified requirements.
- It is the duty of the user to use only faultless tools, for the correct purpose, and protect patients and third parties from all dangers.
- Higher temperature might cause the sense of oral surface tissue unadaptable. Stop using the handpiece if overheated.

4. Assembly and Removal


- **Assembly**
Push the Contra-angle/Straight handpiece onto the motor and turn it until it engages audibly. Check the connection to the motor is secure.
- **Removal**
Remove the Contra-angle / Straight handpiece from the motor by pulling in an axial direction.

Do not assemble / remove during the operation!

5. Forward / Reverse Rotation of Air Motor


- Forward: Turn the FWD/REV Switch Ring to "F" Mode.
- Reverse: Turn the FWD/REV Switch Ring to "R" Mode.
- When FWD/REV Switch Ring is positioned in the middle between "F" and "R", the air motor stops operating.

6. Water and Chip Air Supply


When water is required at the burr, connect the spray tube attachment to the water supply nozzle of the motor and connect the other end to the water supply nozzle of the Straight or Contra-angle Handpiece. Trim the tube to the required length.

* Spray tube attachment for Straight handpiece is an option.


7. Changing Bur

- **Straight Handpiece**
- 1. Turn the Burr Lock Ring anti-clockwise (black arrow) to open the chuck. Insert the burr up to the limit stop.
- 2. Turn the Burr Lock Ring clockwise (white arrow) until it clicks. Check that the handpiece rotary chuck is closed and the burr has engaged.
- 3. Check that the burr is held securely by pulling on it in an axial direction.
- 4. Turn the Burr Lock Ring anti-clockwise (black arrow) to open the chuck. Remove burr from the chuck.

- **Contra-Angle Handpiece**
- 1. Push the latch forward and insert the burr to the center of spindle and turn it until it snaps to place.
- 2. Push the latch forward and take out the burr.

- Use only rotary instruments which are in perfect condition. Follow the operating instructions of the manufacturer.
- Insert rotary instruments only when the Contra-angle / Straight handpiece is stationary.
- Do not interfere with the running or slowing down of rotary instruments.
- Do not activate the rotary chuck of the handpiece or press latch of the Contra-angle during operation or slowing down.

8. Cleaning

- **Contra-angle / Straight handpiece**
 - Remove the rotary instrument.
 - Remove the Contra-angle / Straight handpiece from the motor.
- **Air motor**
 - Remove the transmission instrument from the air motor and set speed control ring to F or R.
 - Remove the air motor from the supply hose .

1. Use clean cloth with mild cleaning solution to clean the handpiece / motor surface carefully, especially groove, and remove the residues.
2. Do not immerse handpiece / motor in disinfectant solution or water. Do not clean by ultrasonic process.
3. Cleaning solution must be non-erosive and harmless to human. Solution with 80% ethyl alcohol is recommended. Please also follow recommendation of the cleaning solution manufacturer for suggested usage.

9. Lubrication

- **Air motor:**
 1. Put Shaft Spray Nozzle(※0PL0044) into nozzle of the lubrication oil can.
 2. Inject the oil into the drive air pipe for 2~3 seconds, at the back end of the motor.
- **Contra-angle / Straight Handpiece:**
 1. Put E-type spray nozzle(※0PL0114) into nozzle of the lubrication oil can.
 2. Connect handpiece with the nozzle and inject lubrication oil for 2~3 seconds.
 3. Clean and lubricate at least once a day after each treatment or before high pressure and high temperature sterilization.

- Use dental handpiece lubrication oil for maintenance, any unspecified oil not intended for lubricating handpieces could shorten the lifetime and harm human body.
- Be sure to lubricate the handpieces prior to sterilization and run the handpieces for 20~30 seconds to emit residual oil and dirt.
- After unpacking and prior to first use, be sure to lubricate with oil and run the handpiece in low rotation speed for 30 seconds and in full speed for 10~20 seconds to emit residual oil.
- Warnings about lubricant hazards: selected spray can labeling and packaging should comply with other Federal requirements for transport and use.

10. Sterilization

- Place handpiece in single-use sterile bag and proceed steam gravity sterilization. Do not open the sterile bag before using in order to avoid contact infection.
- Use qualified sterilization equipment specified by the authority/government for sterilization procedures.
- Do not use chemical steam nor immerse handpiece or motor in disinfectant solution, otherwise could damage the handpiece and cause severe consequences.
- Sterilization time is a minimum of 4 minutes at 134°C, and drying time is 15 minutes (according to EN 13060).

11. Disclaimer


- Once this product is damaged (ie deformed caused by hit) or improperly used (ie sterilized above suggested temperature), be sure not to continue using the handpiece.
- Before operation, always make sure the handpieces and accessories meet operational requirements (ie drive air pressure, water pressure and lubrication, etc)
- Information in this operation manual were carefully reviewed. We cannot assume any liability for any damage arising from improper use of this product or from use of this manual. Please contact our dealer or us with proof of purchase whenever you have problems such as abnormal vibration, exceeding temperature or unusual noise ...etc.
- Products are subject to change without prior notice.
- Products must be operated by qualified practitioners and in dental treatment. Use for other purposes could result in dangers.

CAUTION: U.S. Federal law restricts the use of this device to licensed professionals. (Rx only)

12. Quality and Service

All of our products have passed strict QC procedures before delivery. Please note that improper use or not using the recommended parts and accessories or modifying our products without our approval may decrease the lifetime of the product or cause damage. Please contact our dealers or us with proof of purchase for consultation whenever you have any questions.

Dealer

Manufactured by
TTBIO CORP.
2F., NO., 6TH ROAD INDUSTRY
PARKCHUG. TAIWAN 40755 TEL: 886-4-23595958
FAX: 886-4-23596893
E-mail: medical@ttbio.com
Website: www.ttbio.com

EC REP

Thunder Tiger Europe GmbH
Rudolf-Diesel-Str. 1,
D-86453 Dasing
Germany

* Contents are subject to change without prior notice.

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