


Quiet
Durability
Small Head
Lightweight

VIGOR C36
Contra-angle Handpiece



VIGOR S15
Straight Handpiece



Vigor

C36 / S15
A65 / A68
A65Z

Classic

- GERMAN MADE CERAMIC BALL BEARING
- HIGH STRENGTH GEAR TRANSMISSION
- HIGH OUT-PUT ALLOY CYLINDER AND COMPOSITE PISTON

Vigor SERIES

DURABILITY, SMALL HEAD, AND QUIET

LOW SPEED

RELIABILITY

- 1:1 direct drive transmission ratio
- 25,000 allowable max. speed (rpm)
- Small head design

CA burr spec.

External coolant delivery

Latch chuck type



HP burr spec.

Rotary chuck system for changing the bur easily and quickly



- 1:1 direct drive transmission ratio
- 40,000 allowable max. speed (rpm)
- External coolant delivery
- Low noise and vibration extend the service life of bearing

MODEL

VIGOR C36	C36- 08010143	Contra-angle Handpiece
VIGOR S15	S15- 08010071	Straight Handpiece
VIGOR A65 / A68	A65- 08010072 A68- 08010098	Air Motor
VIGOR A65 Z	A65 Z- 08010185	Air Motor Kit

*Vigor A67 Air Motor without coolant line is available.
E-mail/call us for more product information.

VIGOR A65 / A68 AIR MOTOR

Midwest 4 hole connection
22,000 allowable max. speed (rpm)



- External coolant delivery (A65)



- Internal coolant delivery (A68)

VIGOR A65 Z AIR MOTOR KIT

High power, long life
Smoothly control forward and reverse operating power
Quiet operation, light weight and easy to use



(Vigor A65 + Vigor S15 + Vigor C36) Air Motor Kit

External Type Single Spray
Ample water mists with accurate spray angle.

Durability
Ergonomically Chuck Design.